AMS MONARC®

subfascial hammock

Instructions for Use

AMS MONARC® subfascial hammock Instructions for Use	page 1
AMS MONARC® hamac sous-urétral Mode d'emploi	page 9
AMS MONARC® subfasziale Hängematte Gebrauchsanleitung	Seite 19
AMS MONARC® sostegno sottofasciale Istruzioni per l'uso	pagina 29
AMS MONARC® hamaca subfascial Instrucciones de uso	página 39
AMS MONARC® rede subfascial Instruções de utilização	página 49

European Union (EU) Authorized Representative:

American Medical Systems Europe B.V.

Straatweg 66H
3621 BR Breukelen
The Netherlands
Tel: +31 346 258 100
Fax: +31 346 258 130

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AMS MONARC[®] *subfascial hammock* Instructions For Use

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

These instructions are recommended for general use of this device for the treatment of stress urinary incontinence caused by urethral hypermobility and/or intrinsic sphincter deficiency (ISD). Variations in use may occur in specific procedures due to individual technique and patient anatomy.

DESCRIPTION

The Monarc[®] subfascial hammock is a sterile, single use procedure kit consisting of:

- Two stainless steel, curved needle passers. The tip portion of each needle passer is configured to allow for secure placement of the connectors. Each needle passer has a plastic handle attached.
- An assembly including one piece of loosely knitted polypropylene mesh, two removable plastic insertion sheaths attached to the mesh, and two connectors attached to the insertion sheaths. The polypropylene mesh is constructed of polypropylene monofilament that is pre-cut to 1.1 cm width x 50 cm length. An absorbable tensioning suture is threaded into the length of the mesh to allow for tensioning adjustment of the mesh after placement in the patient is achieved. The two plastic sheaths overlap during placement. The plastic covering also affords convenient travel of the mesh through the tissue. The connectors are attached to the tip portions of the Monarc needle passers during the procedure. The loosely knitted polypropylene mesh is intended to remain in the body as a permanent implant. The mesh component is not absorbed or degraded by the action of tissue in-growth or tissue enzymes.

INDICATIONS

The Monarc subfascial hammock is intended for the placement of a suburethral mesh for the treatment of female stress urinary incontinence (SUI) resulting from urethral hypermobility and/or intrinsic sphincter deficiency.

INSTRUCTIONS FOR USE

Overview

The Monarc procedure can be carried out under local, regional or general anesthesia. The dissection required is minimal. A small incision is made in the anterior vaginal wall followed by periurethral dissection.

Two small stab skin incisions are also made along the medial edge of the obturator foramen at the base of the adductor longus tendon attachment, approximately at the level of the clitoris. Needles are passed from skin to vaginal incisions. The connectors are attached to each needle end; needles are retracted and mesh is positioned without tension.

Procedure

PATIENT PREPARATION


1. Patient should be placed in a modified dorsal lithotomy position with hips flexed, legs elevated in stirrups and buttocks even with edge of table.
2. Ensure the bladder is empty, a catheter is not required during the procedure but may aid in identifying the urethra during the procedure.
3. Place a weighted vaginal retractor or use other suitable vaginal retraction, if desired.

INCISION AND VAGINAL DISSECTION

4. On the anterior wall of the vagina, draw a mark to determine the location of a suburethral incision approximately at the level of the mid-urethra. Variations may occur in specific incisions due to individual technique and patient anatomy.
5. Infiltrate incision site with saline if desired.
6. Place an Allis forceps on the incision margin to expose the incision.
7. Incise the vaginal wall and extend the dissection laterally (periurethral) with a Metzenbaum scissors. Dissect the

periurethral attachment to the vagina.

8. Blunt or sharp dissection may be used. With sharp dissection, insert the tip of the Metzenbaum scissors laterally, spread and advance the scissors until the tip of the scissors touches the interior portion of the inferior pubic ramus (about 1-1.5 cm). Do this bilaterally.
9. The vaginal dissection should be large enough for a fingertip to enter in both directions.
10. Before making skin incisions, palpate the edge of the ischiopubic ramus beginning at the level of the vaginal incision; continue palpating the edge of the bone moving cephalad toward the level of the clitoris and note where the adductor longus tendon inserts into the pubic ramus. Just inferior to this location, and lateral to the bony edge, will be the mark for the skin incisions. Repeat on patient's contralateral side. Confirm that both marks lie in a straight line approximately at the level of the clitoris. Make a small stab incision over these marks.

Small vertical incisions over the obturator foramen

11. The Monarc package can now be opened. Be certain that the package integrity has not been compromised in shipping, and inspect the device for damage.

NEEDLE INSERTION THROUGH THE OBTURATOR FORAMEN

12. Remove Monarc needles from package and grasp the needle designated for the patients left side from the package. Grasp the needle handle with right hand. (See illustration B)

- Point the needle tip perpendicular to the skin with tip in the incision.
 - Place index finger from left hand in the vaginal incision, until fingertip palpates the bone.
13. The thumb from the left hand should be on the outside curve of needle to control the needle movement as it perforates the obturator membrane and muscle.


Illustration A


Illustration B

14. The left thumb pushes the needle through the obturator muscle and membrane. The needle tip penetrates until resistance against the needle stops – about 3-4 cm, depending upon the patient's anatomy.
15. The needle shaft and handle should be positioned at a 45° angle to the patient's vertical axis and close to the patient's body.
16. Rotate the needle handle, moving the needle tip and curve around the posterior surface of the ischial pubic ramus toward the vaginal incision and index finger.
17. If needle tip hits the pubic bone during rotation, retract needle, penetrate needle tip beyond initial insertion depth and rotate toward the vaginal incision.
18. The index finger tip will meet the needle tip as it moves around the pubic ramus. If the needle tip can not be located, then the needle should be withdrawn to just behind the pubic ramus and advanced again.
19. Guide the needle tip with the index finger towards the vaginal incision until the needle tip extends through the incision.

MESH ATTACHMENT AND POSITIONING

20. Attach the connector that is pre-attached to the sheath enclosing the loosely knitted polypropylene mesh to the needle exiting the vaginal incision. The blue markings on

the plastic sheath should be facing outward, away from the urethra. Once the connector is attached to the needle end, reverse rotate the needle back through the skin incision, pulling the connector and associated plastic insertion sheath and mesh into position.

21. Cut the insertion sheath and mesh below the blue mark on the end portion of the plastic sheath and discard the needle, attached connector, sheath end and mesh end. This step allows the sheath to slide freely relative to the mesh.
22. Repeat steps 12-21 on patient's right, or contra-lateral side. (See illustration A)
23. Cystoscopy is not required but can be done at this point at surgeon's discretion.

TENSIONING

24. If vaginal retraction has been used during needle placement, it should be removed when adjusting the tension of the mesh.
25. The mesh is properly tensioned by placing a Metzenbaum scissors or small instrument between the mesh and urethra. To remove the plastic sheaths enclosing the mesh, grasp the ends of the plastic sheaths that have been cut to disconnect the needle and connector. The plastic sheaths are opened and exiting the skin incisions.
26. The mesh is be positioned under the urethra without tension.
27. Once the plastic sheaths have been removed, further adjustment should be minimized.
28. If the mesh needs to be loosened; place an instrument between the mesh and the urethra and pull down, or away from the urethra until proper tension is achieved.
29. If the mesh needs to be tightened, grasp the mesh and tensioning suture exiting the skin incision on each side with a hemostat. Roll around hemostat to improve grip and pull up or out to tighten until proper tension is achieved.
30. Trim the mesh at the level of the subcutaneous tissue and close incisions.
31. Close the vaginal incision.

IMMEDIATE POST-OPERATIVE CARE

32. A catheter and/or vaginal pack can be used at the discretion of the surgeon. Remove prior to discharge.
33. Antibiotic prophylaxis should be given.
34. The ability of the patient to empty the bladder should be confirmed prior to discharge.

CONTRAINDICATIONS

- Do not perform this procedure on pregnant patients.

WARNINGS AND PRECAUTIONS

General

- The risks and benefits of performing the Monarc procedure in women planning future pregnancies should be carefully considered.
- Take special care in cases of bladder prolapse because of anatomical distortion. If the patient requires a cystocele repair, it should be done prior to Monarc procedure. Cystoscopy for these patients can be considered at the doctor's discretion.
- The risks and benefits of performing the Monarc procedure on patients with blood coagulation disorders should be carefully considered.
- The risks and benefits of performing the Monarc procedure on patients with compromised immune systems or any other conditions that would compromise healing should be carefully considered.
- The risks and benefits of performing the Monarc procedure on patients with renal insufficiency and upper urinary tract obstruction should be carefully considered.
- Users should be familiar with surgical procedures and techniques involving non-absorbable meshes before using the Monarc device.
- Users should be familiar with surgical technique for bladder slings and suspensions before employing the Monarc device.
- Users should note the importance of placing the mesh **without tension** under the mid-urethra.
- Good surgical practice should be followed for management

of contaminated or infected wounds.

- Use of the Monarc procedure should only be considered for patients whom the physician determines are adequate surgical risks.
- Vaginal or urinary tract infection should be treated prior to implantation.
- Acute inflammatory tissue reaction and transitory local irritation has been reported with the use of the absorbable suture.

Procedural

- Do not insert the needles with an approach that deploys the needle from the vagina to the obturator foramen.
- It is important to verify the tension and placement of the mesh prior to closure.
- Take care to avoid vessel perforation. Observe patient for any signs of retropubic or periurethral bleeding.
- Take care to avoid perforation of the urethra, bladder or bowel during needle placement.
- Do not remove the plastic sheath until the mesh is in its desired position.
- Do not contact the polypropylene sling with any staples, clips or other instruments which may damage the mesh.

Post Procedure

- If subsequent infection occurs, the entire mesh may have to be removed or revised.
- The patient should be counseled that future pregnancies may negate the effects of the surgical procedure and the patient may again become incontinent.
- Patients should be counseled on abstaining from heavy lifting, exercise and intercourse for a minimum of 4 weeks. Patients can return to other normal daily activities at the physician's discretion, often one to two weeks.
- If dysuria, bleeding or other problems occur, the patient should be instructed to call the surgeon immediately.
- Length of catheter use or vaginal pack should be minimized.

Device Related

- Do not open the secondary compartment of the external pouch. It contains a drying agent (desiccant) that is included as part of the pouch packaging and is not for use with the device. Open the external pouch from the end with the angled, peelable seal only.
- Do not re-sterilize or reuse this device. The Monarc device is intended for single use only. No portion of this procedure kit is reusable.
- Non-functional instruments should not be used and should be returned to AMS.
- Do not use any part of the Monarc procedure kit beyond the indicated expiration date.
- Do not use the Monarc procedure kit if the package is opened or damaged, as sterility may be compromised.
- Store the Monarc procedure kit in a clean, dry, dark area at room temperature.

ADVERSE EVENTS

- As with all implants, local irritation at the wound site and/or a foreign body response may occur.
- Tissue responses to the implant could include vaginal extrusion, erosion through the urethra or other surrounding tissue, migration of the device from the desired location, fistula formation and inflammation. The occurrence of these responses may require removal of the entire mesh.
- Like all foreign bodies, the polypropylene mesh may potentiate an existing infection.
- Excess tension may cause temporary or permanent lower urinary tract obstruction and retention.
- Known risks of surgical procedures for the treatment of incontinence include pain, infection, erosion, device migration, complete failure of the procedure resulting in incontinence and mild to moderate urinary incontinence due to incomplete support or overactive bladder.

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All rights reserved. Printed in USA.
American Medical Systems, Inc.
© 2009 American Medical Systems, Inc.