Dura-II

Concealable Penile Prosthesis

Instructions for Use

Dura-II™ Concealable Penile Prosthesis Instructions for Use	page 1
Prothèse pénienne dissimulable Dura-II™ Mode d'emploi	page 9
Verbergbares Dura-II™ Schwellkörperimplantat Gebrauchsanleitung	Seite 17
Protesi peniena occultabile Dura-II™ Istruzioni per l'uso	pagina 25
Prótesis peniana disimulable Dura-II™ Instrucciones de uso	página 33
Prótese Peniana Ocultável Dura-II™ Instruções de Utilização	página 41

American Medical
Systems

**European Union (EU)
Authorized Representative:**

**American Medical Systems
France S.A.S.**

19 avenue du Québec
91951 Courtaboeuf Cedex
France
Tel 33 1 69 59 97 00
FAX 33 1 69 59 97 29

Contact List:

American Medical Systems U.S.A.

10700 Bren Road West
Minnetonka, MN 55343
U.S.A.
Tel 952 933 4666
FAX 952 930 6157

**American Medical Systems
Australia Pty.**

Suite 2, Level 2
460 Church Street
North Parramatta 2150
Australia
Tel 61 2 8838 8800
FAX 61 2 9890 9899

**American Medical Systems
Benelux B.V.B.A.**

Toekomstraat 10
1930 Zaventem
Belgium
Tel 32 2 720 33 93
FAX 32 2 721 77 71

**American Medical Systems
Canada Inc.**

P.O. Box 461
Guelph, Ontario N1K 1B6
Canada
Tel 519 826 5333
FAX 519 821 1356

**American Medical Systems
Deutschland GmbH**

Kronenstraße 3
D-10117 Berlin
Germany
Tel 49 30 20 64390
FAX 49 30 20 643999

**American Medical Systems
France S.A.S**

19 avenue du Québec
91951 Courtaboeuf Cedex
France
Tel 33 1 69 59 97 00
FAX 33 01 69 59 97 29

**American Medical Systems
Ibérica S.L.**

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DEVICE DESCRIPTION

The Dura-II concealable penile prosthesis consists of a pair of cylinders which are surgically implanted, one into each corpus cavernosum, to provide penile rigidity. Each cylinder consists of three components: a proximal tip, a distal tip, and a body cylinder (see Figure 1). All components consist of implantable medical-grade materials.

The proximal and distal tips used with the device are molded from silicone elastomer and have an integral metal tip pin. Set screws are used to attach these tips to the prosthesis body.


Figure 1. Dura-II Construction

The body cylinder consists of a central section of articulating plastic segments. A cable extends through the central section of these articulating segments and is connected on each end to a spring which is encased in a metal housing. The entire body cylinder is covered with a durable sheath.

The springs apply tension to the cable on either end of the body cylinder to allow sufficient movement of the cable so that the cylinder can be flexed. The articulating segments, under tension by the cable and spring assembly, provide the desired rigidity and have sufficient surface friction to maintain the position in which they are placed. These features allow the implanted Dura-II penile prosthesis to be positioned as desired for concealment or for intercourse.

The Dura-II is available in 10 and 12 mm diameters with an assortment of matching diameter proximal and distal tips.

SIZES OF AVAILABLE COMPONENTS

10 mm Diameter

1 pair body cylinders, 13 cm long, with hex wrench and set screws.

9 pair distal tips

1 cm*
2 cm*
3 cm
4 cm
5 cm
6 cm
7 cm
8 cm
9 cm

9 pair proximal tips

1 cm*
2 cm*
3 cm
4 cm
5 cm
6 cm
7 cm
8 cm
9 cm

*1 and 2 cm tips can be used in either proximal or distal orientation.

12 mm Diameter

1 pair body cylinders, 13 cm long, with hex wrench and set screws.

9 pair distal tips

1 cm*
2 cm*
3 cm
4 cm
5 cm
6 cm
7 cm
8 cm
9 cm

9 pair proximal tips

1 cm*
2 cm*
3 cm
4 cm
5 cm
6 cm
7 cm
8 cm
9 cm

*1 and 2 cm tips can be used in either proximal or distal orientation.

INDICATIONS

The Dura-II concealable penile prosthesis is indicated for the surgical treatment of impotence.

CONTRAINDICATIONS

If the patient's total intracorporal and pubis to mid-glans lengths are outside those shown in the Dura-II concealable penile prosthesis instructions or the separate sizing chart, use of the prosthesis is contraindicated.

Other than the usual contraindications for anesthesia and surgery, no general contraindications for the use of the Dura-II penile prosthesis or any other penile prosthesis are known. However, the patient's age and state of medical or psychological health could affect the successful implantation and use of his prosthesis.

WARNINGS AND RISKS

Penile prosthesis implantation may result in the loss of any natural or latent spontaneous erectile capacity.

Potential complications include, but are not limited to: infection, prolonged postoperative pain or discomfort, tissue sensitivity reactions, loss of penile sensation, migration or erosion of the prosthesis, urinary obstruction or other complications related to prosthesis migration. Mechanical or performance failures which could render the prosthesis inadequate for sexual intercourse or concealment may occur.

These complications may necessitate surgical revision or removal of the prosthesis.

The risk of surgical revision or implant removal is common to all implants. Generally, surgical revision or removal of penile implants are performed to address other complications. However, the patient may elect to have the device removed due to patient dissatisfaction unrelated to an identifiable safety or effectiveness problem. Removal of an implanted prosthesis, for any reason, without timely reimplantation of a new prosthesis may substantially complicate reimplantation or may make it inappropriate.

The patient must be informed a penile implant is subject to wear and is not considered a lifetime prosthesis. The patient should be made fully aware that additional surgery for replacement and removal may be necessary.

The Dura-II prosthesis is designed to be implanted as a pair of matched cylinders. A single implanted cylinder may not be adequate to achieve sexual intercourse and may have an adverse effect on the reliability of the device.

PRECAUTIONS

Direct applications of surgical instruments to the prosthesis may result in damage, rendering the device unsuitable for implantation.

Repeated excessive bending of Dura-II body cylinders during insertion through a peno-scrotal incision may damage the prosthesis and shorten its usual life.

Distal tip trimming is never recommended.

Proximal tips may be trimmed up to 1 cm.

Do not sand any trimmed edges of the prosthesis. The sanding grit cannot be adequately removed in a sterile environment. Residue of sanding grit on the prosthesis can precipitate a severe in vivo inflammatory response.

Implantation of a penile prosthesis which has been in previous contact with, or contaminated by body tissue or fluid, regardless of intervening cleaning or sterilization is prohibited.

Careful intraoperative sizing is required to ensure device operation and to minimize the occurrence of sizing-related complications.

When performing the sub-coronal approach on uncircumcised patients, surgeons may wish to perform a circumcision to reduce the risks of post-operative complications associated with this approach.

Pre-operative physician consultation and evaluation is recommended for any penile prosthesis implant candidate.

Adequate patient education is required to ensure safe, effective patient use.

The Dura-II penile implant is a mechanical device which incorporates articulating segments. Once characteristic of articulating segments is vibration during movement. In the Dura-II this may be detected either physically or audibly by some patients.

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.

The Dura-II concealable penile implant is provided sterile for single use only. **DO NOT RESTERILIZE.**

INSTRUCTIONS FOR USE

Instruments

Assemble the typical surgical instruments required to perform a sub-coronal or peno-scrotal penile prosthesis implantation. In addition, a penile prosthesis sizer set, such as the Penile Prosthesis Sizer Set, and a set of corporal dilators are required. An eyelid retractor may be useful with the sub-coronal approach.

Implantation Protocol

1. **Pubis to Mid-Glans Length Measurement** – Prior to making any incisions, stretch the penis to approximate erection. Externally measure the stretched penis in centimeters from pubis to mid-glans by lightly resting the sizer on the pubic area. *Do not press the sizer into the skin at the base of the penis.* Record this measurement on the implant registration form supplied with the prosthesis body.
2. **Surgical Incision** – A sub-coronal surgical approach is recommended. With this approach, perform standard corporotomies. When a peno-scrotal approach is used, corporotomies 5 to 7 cm must be made to allow for prosthesis insertion. If difficulties are encountered, corporotomies should be extended to facilitate insertion of the prosthesis. When the sub-coronal approach is used, the corporotomies must be long enough to allow the distal tips of the prosthesis to be inserted under the glans.
3. **Corporal Dilation/Prosthesis Diameter Selection** – Dilate both corpora proximally and distally to accommodate prosthesis diameter (10 or 12 mm). Simultaneously insert two same-diameter dilators to verify that the penis can accommodate two similarly sized body cylinders.
4. **Total Intracorporal Length Measurement** – Stretch the penis to approximate an erection. Using the penile prosthesis sizer, measure the stretched corporal length in centimeters in both proximal and distal directions to determine total intracorporal length. Length must be measured for each corpus. Record these measurements on the implant registration form supplied with each pair of body cylinders.
5. **Selection of Prosthesis Components** – Using the measurements for pubis to mid-glans length and total intracorporal length, refer to the sizing chart (see Figure 2) or Dura-II penile prosthesis sizing chart shipped with each kit to select the appropriate proximal and distal tip lengths used.

Sizing Chart for the Dura-II Concealable Penile Prosthesis

Example: Pubis to mid-glans length = 11 cm
 Total intra-corporal length = 20 cm
 Distal tip length = 3 cm
 Proximal tip length = 4 cm

		Pubis to Mid-Glans Length (cm)												
		7	8	9	10	11	12	13	14	15	16	17	18	
Total Intra-corporal Length (cm)	15	1.1	1.1	1.1	1.1									D,P
	16	1.2	1.2	1.2	2.1	2.1								D,P
	17	1.3	1.3	1.3	2.2	3.1	3.1							D,P
	18	1.4	1.4	1.4	2.3	3.2	4.1	4.1						D,P
	19	1.5	1.5	1.5	2.4	3.3	4.2	5.1	5.1					D,P
	20	1.6	1.6	1.6	2.5	3.4	4.3	5.2	6.1	6.1				D,P
	21	1.7	1.7	1.7	2.6	3.5	4.4	5.3	6.2	7.1	7.1			D,P
	22	1.8	1.8	1.8	2.7	3.6	4.5	5.4	6.3	7.2	8.1	8.1		D,P
	23	1.9	1.9	1.9	2.8	3.7	4.6	5.5	6.4	7.3	8.2	9.1	9.1	D,P
	24		2.9	2.9	2.9	3.8	4.7	5.6	6.5	7.4	8.3	9.2	9.2	D,P
	25			3.9	3.9	3.9	4.8	5.7	6.6	7.5	8.4	9.3	9.3	D,P
	26				4.9	4.9	4.9	5.8	6.7	7.6	8.5	9.4	9.4	D,P
	27					5.9	5.9	5.9	6.8	7.7	8.6	9.5	9.6	D,P
	D,P	D,P	D,P	D,P	D,P	D,P	D,P	D,P	D,P	D,P	D,P	D,P	D,P	

Figure 2. Sizing Chart

Note: If the total intracorporal measurements vary from the left and right corpus, redilate both corpora and remeasure. If the measurements still vary, select the shorter length implant.

Key: First number is distal tip length.
 Second number is proximal tip length.

6. **Assemble Prosthesis** – The body cylinder is bi-directional. Proximal or distal tips may be applied to either end. To assemble the tips to the body cylinders, use the hex wrench supplied with each pair of body cylinders to **tighten screws securely** (see Figure 3). Insert tip pin fully into receptacle on end of body cylinder. Hold tip in place while tightening set screw with hex wrench. When set screw is firmly seated, continue turning the hex wrench until the wrench shaft slightly flexes.


Figure 3. Tip Application to the Dura-II Prosthesis

7. **Insert Prosthesis** – When using the sub-coronal approach, first insert the proximal end of the prosthesis into the proximal corpora. Then insert the distal tip of the prosthesis under the glans by lifting and retracting the corpus distally with an eyelid retractor. At the same time, press the distal tip of the prosthesis into the corpus under the glans.

When the peno-scrotal approach is used, corporotomies must be of sufficient lengths to allow the prosthesis cylinders to be inserted to minimize use of excessive force and severe bending.

Note: *If the prosthesis does not fit (either too long or too short) once inserted, remove, remeasure pubis to mid-glans and total intracorporal lengths and refer to sizing chart for correct tip lengths.*

If difficulty is encountered while inserting either of the cylinders, the length of the corporotomy may need to be extended.

Close the incision using acceptable surgical technique.

POSITIONING THE IMPLANTED PROSTHESIS

Stabilization in the desired position (downward concealed or upward erect) is accomplished by underlaying or overlaying a finger at the dorsal or ventral base of the penis and curving the penis into the desired position.

IMPLANT REGISTRATION

American Medical Systems, Inc. (AMS) requires implantation of penile prosthesis to be registered as soon as possible after the implant date. In the United States, this is done by completing and returning the AMS Patient Information Form (PIF) supplied with each pair of Dura-II body cylinders. Outside the United States, return the appropriate information to the distributor who sold the prosthesis.

The purpose of the registration of implanted prosthesis is to record a complaint, or to identify and protect the patient in the event of a device recall or warranty claim.

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All rights reserved. Printed in USA.
© 2003 American Medical Systems, Inc.