Straight·In™

Fixation System

Instructions for Use

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**European Union (EU) Authorized
Representative:**

American Medical Systems

Europe B.V.
Straatweg 66H
3621 BR Breukelen
The Netherlands
Phone: 31 346 258 100
Fax: 31 346 258 130

Contact List:

American Medical Systems U.S.A.

10700 Bren Road West
Minnetonka, MN 55343
U.S.A.
Tel: 952 930 6000
Fax: 952 930 6157

American Medical Systems

Australia Pty Ltd.

Unit 39, Building F
16 Mars Road
Lane Cove NSW 2066
Australia
Tel: 61 2 9425 6800
Fax: 61 2 9427 6296

American Medical Systems Canada Inc.

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American Medical Systems Ibérica S.L.

c/Joaquin Turina,
2 Planta primera - Oficina 6
28224 Pozuelo de Alarcón (Madrid)
Spain
Tel: 34 91 799 49 70
Fax: 34 91 715 75 26

American Medical Systems U.K. Ltd.

Capital Court
Capital Interchange Way
Brentford
TW8 OEX
United Kingdom
Tel: 44 (0) 20 8996 3100
Fax: 44 (0) 20 8995 3720

Instructions For Use

STRAIGHT·IN™ Fixation System

DESCRIPTION

The Straight-In™ Fixation System is for soft tissue fixation to bones in the pelvic region by means of bone screws threaded with suture. The disposable, battery operated Straight-In Inserter (Figure 1) is available with a 28 cm shaft and is suitable for both open and laparoscopic approaches. The Inserter may be used with either the AMS Bone Screw with polypropylene suture (Figure 2a) or the AMS Loop-Screw with braided polyester suture (Figure 2b). Both have the same miniature shape and sharp tip to facilitate their penetration into the bone without the need for any pre-drilling. The AMS Loop-Screw has a short polypropylene No. 1 loop through which braided polyester No. 0 suture is threaded. The AMS Bone Screw with polypropylene suture has two pre-attached strands of polypropylene No. 1 suture. The specifications of the bone and loop screws are detailed in Table 1.


Figure 1: The assembled Straight-In Bone Screw Inserter


Figure 2a:
The AMS Bone Screw with Polypropylene Suture


Figure 2b:
The AMS Loop-Screw with Braided Polyester Suture

Table 1: Specifications of the AMS Bone and Loop Screws	
Thread Maximal Diameter:	3.2 mm
Length:	5.5 mm
Material:	Titanium

INDICATIONS

The Straight-In Fixation System is intended for soft tissue fixation to bones in the pelvic region (e.g., sacral) by means of bone screws threaded with suture. It is indicated for use during open or laparoscopic surgical procedures where soft tissue fixation to bone in the pelvic region is needed.

CONTRAINDICATIONS

1. Bone deformation or other pathological conditions of bone, such as severe osteoporosis, which would impair secure bone screw fixation.
2. Compromised immune system.
3. Patients having urinary tract infections or previous or active osteomyelitis.
4. Renal insufficiency and upper urinary tract relative obstruction.

This list is not intended to be an exhaustive list and does not cover the types of risks that are attendant to any open surgical procedure.

PRECAUTIONS

1. The Straight-In Fixation System should be used only by physicians who have been trained on its use.
2. The sterile packaging of the Inserter and AMS bone screws should be inspected for visible damage prior to use. Do not use if damage is suspected.
3. The Straight-In Inserter is a single-patient use device. The screws and suture are for single use only. Neither re-sterilization nor reuse should be attempted.
4. The patient should be counseled to avoid physical strain and lifting for 2-3 months post-procedure.

STERILITY

The Straight-In Inserter, 28 cm shaft, and AMS bone and loop screws with attached suture are supplied sterile (EtO) for single-patient use only. All sterile packaging should be inspected for damage and expiration prior to use. Do not use if it is suspected that sterility has been compromised. Re-sterilization should not be attempted.

COMPLICATIONS

Osteitis pubis and osteomyelitis are rare but potential complications. Known risks associated with current open and laparoscopic surgical procedures for soft tissue fixation to one of the pelvic bones include surgical infection and pain. Other potential complications are procedure associated.

STRAIGHT-IN SHAFT ASSEMBLY PROCEDURE

Remove both the shaft and handle from their packaging. Align the shaft within its outlet, press the Shaft Locking Switch and press the shaft to click it into place. For disassembling the shaft, press the Shaft Locking Switch and remove the shaft.

AMS LOOP-SCREW WITH BRAIDED POLYESTER SUTURE LOADING PROCEDURE

The Inserter Safety Switch should be in the READY position, when passing sutures through the Inserter.

Note: The 28 cm shaft should be attached to the Inserter prior to starting this procedure. The braided suture is covered by a stiff plastic tube (Figure 2b) to facilitate threading through the shaft.

1. Move the safety switch of the Inserter to the READY position.
2. Thread the suture/tubing (end opposite the screw) through the shaft until it protrudes from the suture exit hole (Figure 3).
3. Gently pull on the protruding end of the suture/tubing until the screw approaches the screw outlet.
4. Holding the screw by its holder, insert the screw base into the outlet, ensuring that the screw base and outlet are in proper alignment.
5. Press the screw base into the outlet until the screw clicks into place. The screw holder will dislodge and may be discarded.
6. Pull the tubing to slide it off the suture and out the exit hole.
7. Pull on the protruding ends of the suture and move the safety switch to the SAFE position.
8. The Straight-In is now ready for use.
9. This procedure is repeated for each screw loading.


Figure 3

AMS BONE SCREW WITH POLYPROPYLENE SUTURE LOADING PROCEDURE

The Inserter Safety Switch should be in the READY position, when passing sutures through the Inserter.

Note: The shaft should be attached to the Inserter prior to starting this procedure.

1. Move the safety switch to the READY position.
2. Thread the suture (end opposite the screw) through the shaft until it protrudes from the suture exit hole (Fig. 3).
3. Gently pull the protruding end of the suture until the screw approaches the screw outlet.
4. Holding the screw by its holder, insert the screw base into the outlet, ensuring that the screw base and outlet are in proper alignment.
5. Press the screw base into the outlet until the screw clicks into place. The screw holder will dislodge and may be discarded.
6. Pull on the protruding ends of the suture and move the safety switch to the SAFE position.
7. The Straight-In is now ready for use.
8. This procedure is repeated for each screw loading. The suture ends should not be cut until confirmation of screw insertion into bone is determined.

SCREW INSERTION PROCEDURE

1. Place the tip of the screw perpendicular to the bone.
2. Move the safety switch to the READY position.
3. Press the Inserter's operating button until the screw is completely inserted. Then release the operating button and remove the Inserter.
4. The lot numbers of the screws used in the procedure should be recorded on the patient's chart (pre-printed peel-off labels are provided on the screw package).
5. **For U.S.A. procedures only:** At the end of the procedure, complete the Patient Information Form provided with each Inserter package and send it to American Medical Systems in the enclosed envelope.
6. This procedure is repeated for each screw insertion.

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

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American Medical Systems, Inc.