AMS BioArc® TO Subfascial Hammock with Pre-Connected InteXën® LP (AMS Collagen Dermal Matrix) Instructions for Use

AMS BioArc® TO page 1
Subfascial Hammock with
Pre-Connected InteXën® LP
(AMS Collagen Dermal Matrix)
Instructions for Use

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American Medical Systems Australia Pty Ltd.

Unit 39, Building F
16 Mars Road
Lane Cove NSW 2066
Australia
Tel: +61 2 9425 6800
Fax: +61 2 9427 6296

American Medical Systems Canada Inc.

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American Medical Systems Ibérica S.L.

c/Joaquín Turina, 2 Planta primera – Oficina 6
28224 Pozuelo de Alarcón (Madrid)
Spain
Tel: +34 (0) 91 799 49 70
Fax: +34 (0) 91 715 75 26

American Medical Systems U.K. Ltd.

Capital Court
Capital Interchange Way
Brentford
TW8 0EX
United Kingdom
Tel: +44 (0) 20 8996 3100
Fax: +44 (0) 20 8995 3720

AMS BioArc® TO Subfascial Hammock with Pre-Connected InteXēn® LP (AMS Collagen Dermal Matrix)


Instructions For Use

These instructions are recommended for general use of this device for the treatment of stress urinary incontinence caused by urethral hypermobility and/or intrinsic sphincter deficiency (ISD). Variations in use may occur in specific procedures due to individual technique and patient anatomy.

DESCRIPTION

The BioArc TO Subfascial Hammock with Pre-Connected InteXēn LP is a sterile, single use procedure kit consisting of:

- Two stainless steel, curved, needle passers. The tip portion of each needle passer is configured to allow for secure placement of the connectors. Each needle passer has a plastic handle attached.
- An assembly including one piece of loosely knitted AMS polypropylene mesh, with attached 1.8 x 7 cm (after hydration) collagen dermal graft, two removable plastic insertion sheaths attached to the mesh, and two connectors attached to the sheath. The graft is sterile, non-perforated processed porcine collagen dermal matrix that is lyophilized.

The AMS mesh is constructed of polypropylene monofilament that is precut to 1.1 cm width and includes a 1.8 x 7 cm (after hydration) collagen dermal graft in the center, to be placed suburethrally. A polypropylene tensioning suture is threaded into the length of the polypropylene mesh to allow tightening of the subfascial hammock after placement in the patient. A polyester loosening suture with ID tab is also threaded through the mesh to allow loosening of the subfascial hammock after placement in the patient without distorting the graft. The plastic insertion sheaths enclose the mesh and are designed to protect the mesh during placement. The connectors are attached to the tip portion of the BioArc TO needle passers during the procedure. The loosely knitted AMS polypropylene mesh and collagen dermal graft is intended to remain in the body as a permanent implant. The mesh component is not absorbed or degraded by the action of tissue ingrowth or tissue enzymes.

INDICATIONS

The BioArc TO subfascial hammock with Pre-Connected AMS Collagen Dermal Matrix is intended for the placement of a suburethral graft for the treatment of female stress urinary incontinence (SUI) resulting from urethral hypermobility and/or intrinsic sphincter deficiency.

Overview

The BioArc TO procedure can be carried out under local, regional or general anesthesia. The dissection required is minimal. A small incision is made in the anterior vaginal wall followed by periurethral dissection.

Two small stab skin incisions are also made along the medial edge of the obturator foramen at the base of the adductor longus tendon attachment, approximately at the level of the clitoris. Needles are passed from skin to vaginal incisions. The connectors are attached to each needle end; needles are retracted and the graft is positioned without tension.

Procedure

1. PATIENT PREPARATION

- 1.1 Patient should be placed in a modified dorsal lithotomy position with hips flexed, legs elevated in stirrups and buttocks even with edge of table.
- 1.2 Ensure the bladder is empty. A catheter is not required during the procedure but may aid in identifying the urethra during the procedure.
- 1.3 Place a weighted vaginal retractor or use other suitable vaginal retraction, if desired.


2. MESH PREPARATION AND GRAFT HYDRATION PROCEDURE

- 2.1 Peel open the outer pouch and aseptically remove the inner pouch.
- 2.2 Only the inner pouch may be placed in a sterile field.
- 2.3 Aseptically remove the sling from the inner sterile pouch.
- 2.4 Fill sterile dish with sufficient amount of rehydration fluid, such as normal saline to submerge the sling with attached collagen dermal graft material. Submerge the graft completely and allow it to soak for a minimum of 5 minutes, but not more than 1 hour.
- 2.5 The graft is rehydrated once it's soft and pliable throughout. The graft may require a gentle massage to complete the rehydration process.

3. INCISION AND VAGINAL DISSECTION

- 3.1 On the anterior wall of the vagina, draw a mark to determine the location of a suburethral incision approximately at the level of the mid-urethra. Variations may occur in specific incisions due to individual technique and patient anatomy.
- 3.2 Infiltrate incision site with saline if desired.
- 3.3 Place an Allis forceps on the incision margin to expose the incision.
- 3.4 Incise the vaginal wall and extend the dissection laterally (periurethral) with a Metzenbaum scissors. Dissect the periurethral attachment to the vagina.

- 3.5 Blunt or sharp dissection may be used. With sharp dissection, insert the tip of the Metzenbaum scissors laterally, spread and advance the scissors until the tip of the scissors touches the interior portion of the inferior pubic ramus (about 1-1.5 cm). Do this bilaterally.
- 3.6 The vaginal dissection should be large enough for a fingertip to enter in both directions.
- 3.7 Before making skin incisions, palpate the edge of the ischiopubic ramus beginning at the level of the vaginal incision; continue palpating the edge of the bone moving cephalad toward the level of the clitoris and note where the adductor longus tendon inserts into the pubic ramus. Just inferior to this location, and lateral to the bony edge, will be the mark for the skin incisions. Repeat on patient's contralateral side. Confirm that both marks lie in a straight line approximately at the level of the clitoris. Make a small stab incision over these marks.
- 3.8 The BioArc TO package with needles can now be opened. Be certain that the package integrity has not been compromised in shipping, and inspect the device for damage.

Small vertical incisions over the obturator foramen

4. NEEDLE INSERTION THROUGH THE OBTURATOR FORAMEN

- 4.1 Ensure the bladder is empty. Remove BioArc TO needles from package and grasp the needle designated for the patient's left side from the package. Grasp the needle handle with right hand. (See illustration B)
 - Point the needle tip perpendicular to the skin with tip in the incision.
 - Place index finger from left hand in the vaginal incision, until fingertip palpates the bone.
- 4.2 The thumb from the left hand should be on the outside curve of needle to control the needle movement as it perforates the obturator membrane and muscle.


Illustration A



Illustration B

- 4.3 The left thumb pushes the needle through the obturator muscle and membrane. The needle tip penetrates until resistance against the needle stops – about 3-4 cm, depending upon the patient’s anatomy.
- 4.4 The needle shaft and handle should be positioned at a 35-45° angle to the patient’s vertical axis and close to the patient’s body.
- 4.5 Rotate the needle handle, moving the needle tip and curve around the posterior surface of the ischial pubic ramus toward the vaginal incision and index finger.
- 4.6 If needle tip hits the pubic bone during rotation, retract needle, penetrate needle tip beyond initial insertion depth and rotate toward the vaginal incision.
- 4.7 The index finger tip will meet the needle tip as it moves around the pubic ramus. If the needle tip can not be located, then the needle should be withdrawn to just behind the pubic ramus and advanced again.
- 4.8 Guide the needle tip with the index finger towards the vaginal incision until the needle tip extends through the incision.

5. SUBFASCIAL HAMMOCK ATTACHMENT AND POSITIONING

- 5.1 Attach the connector, that is pre-attached to the sheath enclosing the polypropylene mesh, to the needle exiting the vaginal incision. The blue markings on the plastic sheath should be facing outward, away from the urethra. Be sure that the graft lies flat and that the mesh is not twisted prior to attaching the second connector. **The connectors cannot be removed once they are snapped into place.** Once the connector is attached to the needle end, reverse rotate the needle back through the skin incision, pulling the connector and associated plastic insertion sheath and mesh into position.

- 5.2 Cut the insertion sheath and mesh below the blue mark on the end portion of the plastic sheath and discard the needle, attached connector, sheath end and mesh end. This step allows the sheath to slide freely relative to the mesh.
- 5.3 Repeat steps 17-26 on patient’s right, or contra-lateral side. (See illustration A)
- 5.4 Cystoscopy is not required but can be done at this point at surgeon’s discretion.
- 5.5 Position the graft under the midurethra without tension. Prior to removing the protective sheaths, make sure that the graft material is centered under the urethra with the suture knots facing the surgeon and away from the urethra.

5.6 After desired placement is achieved, remove both plastic sheaths enclosing the mesh by pulling up on each sheath, one side at a time. To avoid overtightening the subfascial hammock, keep a forceps or other instrument between the graft material and urethra during removal.

5.7 The tensioning suture allows intraoperative adjustment of the subfascial hammock after the plastic sheaths have been removed.


To tighten the subfascial hammock:

Grasp the mesh and tensioning suture exiting the skin incision on each side with a hemostat. The mesh may be rolled around the clamp to improve the grip. Pull up to tighten the subfascial hammock as desired.

To loosen the subfascial hammock:

Use a hemostat or a clamp to pull from each of the hanging loosening sutures. Use the clamps to pull down and loosen the subfascial hammock as desired.

Caution: Do not pull ID tab on loosening suture to loosen subfascial hammock.

6. CLOSURE AND POST-OPERATIVE CARE

- 6.1 Cut one end of each loosening suture and pull ID tab until the entire loosening suture is removed.
- 6.2 Trim the mesh to size at the skin incisions.
- 6.3 Close the skin and vaginal incisions.
- 6.4 A catheter and/or vaginal pack can be used at the discretion of the surgeon. Remove prior to discharge.
- 6.5 Antibiotic prophylaxis should be given.
- 6.6 The ability of the patient to empty the bladder should be confirmed prior to discharge.

CONTRAINDICATIONS

- Do not perform this procedure on pregnant patients.
- Preconnected porcine dermis should not be used on patients with known sensitivity to porcine products or on patients with a history of multiple or serum allergies.

WARNINGS AND PRECAUTIONS

General

- The risks and benefits of performing the BioArc TO procedure in women planning future pregnancies should be carefully considered.

- Take special care in cases of bladder prolapse because of anatomical distortion. If the patient requires a cystocele repair, it should be done prior to the BioArc TO procedure. Cystoscopy for these patients can be considered at the doctor's discretion.
- The risks and benefits of performing the BioArc TO procedure on patients with blood coagulation disorders should be carefully considered.
- The risks and benefits of performing the BioArc TO procedure on patients with compromised immune systems or any other conditions that would compromise healing should be carefully considered.
- The risks and benefits of performing the BioArc TO procedure on patients with renal insufficiency and upper urinary tract obstruction should be carefully considered.
- Users should be familiar with surgical procedures and techniques involving non-absorbable meshes before using the BioArc TO Subfascial Hammock.
- Users should be familiar with surgical technique for bladder slings and suspensions before employing the BioArc TO Subfascial Hammock.
- Users should note the importance of placing the subfascial hammock without tension under the mid-urethra.
- Good surgical practice should be followed for management of contaminated or infected wounds.
- Use of the BioArc TO procedure should only be considered for patients whom the physician determines are clinically acceptable surgical risks.
- Vaginal or urinary tract infection should be treated prior to implantation.

Procedural

- A transverse vaginal incision is not recommended.
- Do not insert the needles with an approach that deploys the needle from the vagina to the obturator foramen.
- It is important to verify the tension and placement of the subfascial hammock prior to closure.
- **Take care to avoid vessel perforation. Observe patient for any signs of retropubic or periurethral bleeding.**
- Take care to avoid perforation of the urethra, bladder or bowel during needle placement.
- Do not remove the plastic sheath until the subfascial hammock is in its desired position.
- Do not contact the polypropylene mesh or graft with any staples, clips or other instruments which may damage the subfascial hammock.

Post Procedure

- If subsequent infection occurs, the entire subfascial hammock may have to be removed or revised.
- The patient should be counseled that future pregnancies may negate the effects of the surgical procedure and the patient may again become incontinent.
- Patients should be counseled on abstaining from heavy lifting, exercise and intercourse for a minimum of 4 weeks. Patients can return to other normal daily activities at the physician's discretion, often one to two weeks.
- If dysuria, bleeding or other problems occur, the patient should be instructed to call the surgeon immediately.
- Length of catheter use or vaginal pack should be minimized.

Device Related

- Do not re-sterilize or reuse this device. The BioArc TO device is intended for single use only. No portion of this procedure kit is reusable.
- Inspect each component of the system prior to use. Non-functional components should not be used and should be returned to AMS.
- Do not use any part of the BioArc TO procedure kit beyond the indicated expiration date.
- Do not implant the AMS polypropylene mesh with collagen dermal matrix without first hydrating the graft.
- Do not use the BioArc TO procedure kit if the package is opened or damaged, as sterility may be compromised.
- STORE AT ROOM TEMPERATURE. DO NOT USE PRODUCT IF TEMPERATURE EXCEEDS 99° F (37° C) OR INDICATOR IS BLACK. Immediately return unusable product to American Medical Systems.

ADVERSE EVENTS

- As with all implants, local irritation and/or dehiscence at the wound site and/or a foreign body response may occur.
- Tissue responses to the implant could include vaginal extrusion, erosion through the urethra or other surrounding tissue, migration of the device from the desired location, fistula formation and inflammation. The occurrence of these responses may require removal of the entire subfascial hammock
- Like all foreign bodies, the polypropylene mesh with collagen dermal graft may potentiate an existing infection.
- Excess tension may cause temporary or permanent lower urinary tract obstruction and retention.
- Known risks of surgical procedures for the treatment of incontinence include pain, infection, erosion, device migration, hematoma, blood transfusion, complete failure of the procedure resulting in incontinence and mild to moderate urinary incontinence due to incomplete support or overactive bladder.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

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All rights reserved. Printed in USA.
American Medical Systems, Inc.
© 2010 American Medical Systems, Inc.