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CAUTION-Federal (USA) law restricts this device to sale by or on the order of a physician. Physicians should be trained in the minimally-invasive placement of surgical mesh devices for treatment of pelvic floor disorders and in management of complications resulting from the procedures.

Training on the use of the AMS Elevate® System is recommended and available from the manufacturer. Please contact your AMS Sales Representative to arrange for this training. Please contact AMS Customer Service for additional instructional materials.

The safety and effectiveness of synthetic mesh or film support in transvaginal surgical procedures to treat pelvic organ prolapse have not been evaluated in a prospective, randomized clinical study.

European Union (EU) Authorized Representative:

American Medical Systems

Europe B.V.

Straatweg 66H
3621 BR Breukelen
The Netherlands
Tel: +31 (0) 346 258 100
Fax: +31 (0) 346 258 130

Contact List:

American Medical Systems U.S.A.

10700 Bren Road West
Minnetonka, MN 55343
U.S.A.
Tel: +952 930 6000
Fax: +952 930 6157

American Medical Systems Australia Pty Ltd.

Unit 39, Building F
16 Mars Road
Lane Cove NSW 2066
Australia
Tel: +61 2 9425 6800
Fax: +61 2 9427 6296

American Medical Systems Canada Inc.

P.O. Box 461
Guelph, Ontario N1H 6K9
Canada
Tel: +519 826 5333
Fax: +519 821 1356

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American Medical Systems Ibérica S.L.

c/Joaquín Turina, 2 Planta primera – Oficina 6
28224 Pozuelo de Alarcón (Madrid)
Spain
Tel: +34 (0) 91 799 49 70
Fax: +34 (0) 91 715 75 26

American Medical Systems U.K. Ltd.

Capital Court
Capital Interchange Way
Brentford
TW8 0EX
United Kingdom
Tel: +44 (0) 20 8996 3100
Fax: +44 (0) 20 8995 3720

AMS Elevate® Prolapse Repair Systems

*Apical and Posterior Prolapse Repair System,
Anterior and Apical Prolapse Repair System*

INSTRUCTIONS FOR USE

The AMS Elevate System is intended for use where the connective tissue has ruptured or for implantation to reinforce soft tissues where weakness exists in the urological, gynecological, and gastroenterological anatomy. This includes but is not limited to the following procedures: pubourethral support, including urethral slings, vaginal wall prolapse repairs including anterior and posterior wall repairs, vaginal suspension, reconstruction of the pelvic floor and tissue repair.

CONTRAINDICATIONS

- AMS Elevate System should not be used in infants, children, pregnant women, or women planning future pregnancies as the mesh will not stretch significantly as the patient grows.
- Do not implant in patients with pre-existing conditions that pose an unacceptable surgical risk.
- AMS Elevate System must not be used following planned intra-operative or accidental opening of the gastrointestinal tract. Use in these cases may result in contamination of the mesh, which may lead to infection that may require removal of the mesh.
- AMS Elevate System should not be used in the presence of active or latent infections or cancers of the vagina, cervix, or uterus **or in patients who have received radiation in the area of treatment.**
- If using AMS Elevate System with AMS InteXēn™ LP™ porcine dermis graft, the product should not be used on patients with a known sensitivity to PORCINE products, or on patients with history of multiple or serum allergies.
- Polypropylene mesh should not be placed in contact with bowel or visceral organs, including the urinary bladder.

WARNINGS

For all Pelvic Floor Repair Procedures:

- Only trained physicians should use this product.
- Patients on anticoagulation agents undergoing surgery must have their anticoagulation therapy carefully managed.
- A digital rectal exam should be performed to detect **or avoid** possible rectal perforation.

- Over tensioning may cause temporary or permanent urinary **or defecatory** obstruction.
- The implanted mesh sling must always be separated from the abdominal cavity by peritoneum.

WARNINGS

For Anterior Pelvic Floor Repairs:

- When performing an anterior prolapse repair procedure, **cystoscopy** should be performed to confirm bladder integrity or detect bladder or ureteral perforation.
- Following a bladder neck suspension procedure, the patient should be counseled that future pregnancies may negate the effects of surgical procedure and the patient may again become incontinent.
- Special attention should be given to patients with renal insufficiency or patient with upper urinary tract obstruction.

PRECAUTIONS

- Users should be well informed of surgical procedures and techniques involving pelvic floor repair and synthetic meshes before employing AMS Elevate System.
- As anatomy of individual patients may vary greatly, for each procedure it is important that the intended planes for device advancement and the intended place for mesh placement are planned and known for each individual patient. Employment of imaging methods before and after mesh placement may aid in proper mesh placement and confirm absence of injury to non-target anatomical structures.
- Avoid placing excessive tension on the mesh implant during placement.
- Do not affix AMS Elevate mesh with any staples, clips, or clamps as mechanical damage to the mesh may occur.
- The product should only be used under the prescription of a physician.
- In patients with compromised immune systems or other conditions that would compromise healing the risks and benefits should be carefully weighed.
- Vaginal or urinary tract infection should be treated prior to implantation.
- Prolapse repair may unmask pre-existing incontinence conditions.
- Prophylactic antibiotics can be administered according to the surgeon's usual practice.
- No data is available on the use of AMS Elevate mesh with tissue adhesives and is not recommended.
- Serious adverse tissue responses may require partial or complete removal of the mesh.
- When passing the needles:
 - Take care in cases of bladder prolapse due to anatomical distortion.

- **Take care to avoid vessel perforation or nerve damage. Observe patient for any signs of bleeding.**

- Do not use beyond indicated expiration date.
- Inspect package for damage. If damaged, do not use as sterility may be compromised.
- For devices with AMS InteXën LP porcine dermis graft, STORE AT CONTROLLED ROOM TEMPERATURE. DO NOT USE PRODUCT IF TEMPERATURE EXCEEDS 99° F (37° C) OR INDICATOR IS BLACK. Immediately return unusable product to American Medical Systems.

POTENTIAL ADVERSE EVENTS

- Pain/Discomfort/Irritation
- Inflammation
- Infection
- Mesh erosion
- Mesh extrusion
- Fistula formation
- Foreign body reaction to mesh (includes reaction to porcine products)
- Adhesion formation
- Urinary incontinence
- Urinary retention/obstruction
- Voiding dysfunction
- Recurrent prolapse
- Contracture
- Wound dehiscence
- Nerve damage
- Perforation of vessels, nerves, bladder, ureter, colon, and other pelvic floor structures
- Hematoma
- Dyspareunia

POST-OPERATIVE CARE

- Post-operatively the patient should be advised to refrain from intercourse, heavy lifting and/or exercise (e.g. cycling, jogging) until the physician determines when it is suitable for the patient to return to her normal activities.
- If dysuria, bleeding, dyspareunia, severe pain, defecatory or other problems occur, the patient should be instructed to call the surgeon immediately.

PRODUCT DESCRIPTION

The AMS Elevate is a total Trans-Vaginal System consisting of stainless steel needles used to pass monofilament polypropylene mesh sling for placement into the pelvic floor. Various configurations of the polypropylene mesh sling and needle passers are available depending on the surgical procedure being performed. The mesh sling is also available with AMS InteXën LP porcine dermis graft.

HOW SUPPLIED AND STORAGE

WARNING: Contents supplied STERILE. Do not use if sterile barrier is damaged. If damage is found, call your AMS representative.

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Do not use product beyond indicated expiration date.

Recommended storage conditions: Store product at controlled room temperature.

For devices with AMS InteXën LP porcine dermis graft, STORE AT CONTROLLED ROOM TEMPERATURE. DO NOT USE PRODUCT IF TEMPERATURE EXCEEDS 99° F (37° C) OR INDICATOR IS BLACK. Immediately return unusable product to American Medical Systems.

PRODUCT PREPARATION

- If using AMS Elevate System with AMS InteXën LP, the porcine dermis graft will require re-hydration in sterile saline prior to implantation for a minimum of 5 minutes and no longer than 1 hour. The graft may require a gentle massage to complete the re-hydration process.

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All rights reserved. Printed in USA.
American Medical Systems, Inc.
©2009 American Medical Systems, Inc.